

TITAN STANDARD POWER UNIT

The Titan Standard Power Unit for use with the Titan SGS stapler uses visual and audible status indicators to provide ongoing information to the user and surgical team during the case. Intelligent tissue-sensing technology and enhanced stapling safety features constantly monitor tissue loads and prevent unsafe conditions enabling the surgeon to achieve more consistent results in gastric tissue.




A COMPLETE SOLUTION

Designed for use with the Titan SGS, the stapler utilizes the Standard Bougie™ 38 FR and the bladeless technology of the Standard Trocar™.




TITAN SGS™ EXCLUSIVELY FOR BARIATRIC SURGEONS



ALIGN

The entire stomach tissue is placed into the Titan SGS allowing visualization of the full staple cut line in a single plane.



DIVIDE

Titan SGS fires a single cartridge, completing the stapling in approximately 55 seconds.



ACHIEVE

The ideal tubular surgical sleeve anatomy is a consistent shape that is free of kinks, twists, or spirals.

TITAN SGS SPECIFICATIONS



Description	Titan SGS 23cm powered endoscopic nonarticulating gastric surgical stapler
Order Code	SGS23R
Jaw Length	23cm
Staple Rows ¹	6
Closed Staple Height	1.2 – 2.2mm
EA/BX	3



TITAN STANDARD POWER UNIT SPECIFICATIONS



Description	The Titan Standard Power Unit provides power to the Titan SGS. The SPU01 uses visual and audible status indicators to provide feedback to the user during the case. Reusable (5) years of useful life.
Order Code	SPU01
EA/BX	1

For more information or to place an order: orders@standardbariatrics.com

Device Description: The Titan SGS is a sterile, single-patient-use instrument used for resecting and stapling the stomach during laproscopic and open surgical procedures. **Disclaimer:** Before using any medical device review all relevant package inserts with particular attention to the indications, contraindications, warnings and precautions, and steps for use of the device.

¹Data on file



©2021 Standard Bariatrics, Inc. 4362 Glendale Milford Road, Cincinnati, OH 45242